



Please print at  
100% for the most  
accurate sizing.

**Baxa  
Adapta-Cap™  
Template**




Use as illustrated to determine  
appropriate Adapta-Cap Bottle  
Adapter size.



The smallest opening that the  
screw threads will pass through  
indicates correct Adapta-Cap  
Bottle Adapter size.

**SAVE THIS FOR FUTURE USE.**

Your small friend, indeed®  
**Health Care Logistics**  
INC.



Call Free: 1-800-848-1633 • [www.HealthCareLogistics.com](http://www.HealthCareLogistics.com)

Adapta-Cap is a trademark of Baxa Corporation  
5300-0891 Rev. B

A B C D\* or M\* E F G H I K