

Cardinal Health™ ESTEEM™ Sterile Nitrile Exam Gloves

Delivering protection, flexibility and durability



Powder-free sterile nitrile exam gloves

Count on Cardinal Health for quality exam gloves to help your hands extend patient care.

Features and benefits:

- Extended cuff with a minimum length of 10.6 in.
- Sterile, individually wrapped gloves
- Helps healthcare facilities comply with USP <800> and USP <797> standards for the safe handling of hazardous drugs during sterile compounding procedures¹
- Tested for use with chemotherapy drugs per ASTM D6978^{2,3}
- Beaded cuff assists with donning and helps reduce roll-down
- Textured surface helps improve wet and dry grip

Pairs

Catalog No.	Size	Pkg.
88SNP02S	S	50 pairs/box; 4 boxes/case
88SNP03M	M	50 pairs/box; 4 boxes/case
88SNP04L	L	50 pairs/box; 4 boxes/case
88SNP05XL	XL	50 pairs/box; 4 boxes/case

¹These products help healthcare facilities comply with the USP <797> and USP <800> standard for safe handling of hazardous drugs when used as directed in USP <797> and USP <800>.

²This glove has been tested for resistance to permeation of various chemotherapy drugs per ASTM D6978 "Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs." When chemo drugs are present, gloves selection should be based on the specific types of drugs used. Note: Review safety data sheets (SDS) for the chemicals and drugs used to determine the required level of protection.

³**Warning:** Do not use with Carmustine (3.3 mg/mL) and Thiotepea (10 mg/mL).

- Sterile
- Powder-free
- Not made with natural rubber latex
- Ambidextrous
- Color: Blue
- Minimum fingertip thickness: 4.7 mil (0.12 mm)



ISO-MED

**Product in-Stock
Order while supplies last!**

For more information or to order, contact your ISO-MED representative, call **800-797-1405** or visit **iso-med.com**

